



UNITY™ VALUE
ADDED SERVICE

Installation Services For Laboratory Equipment

Let our service professionals handle your installation. Unity Lab Services will provide you with:

- **Improved Productivity**

Accurate installation scheduled at your convenience, optimizing both equipment and user performance.

- **Peace of Mind**

Your equipment will be installed and configured to our manufacturing specifications ensuring optimal performance while users are trained on equipment operation and routine care.



Installation by the product expert includes:

1. System installation
2. Configuration
3. Functionality test
4. User orientation
5. Customer acceptance



Installation by product experts allows you to focus on what's really important – *your next discovery.*

Unity Lab Services

Installation Services for Laboratory Equipment

Thermo Scientific™ product names and product families covered by these services include:

PRODUCT NAMES

Barnstead™	Lab-Line™
Fisher Isotemp®	Lindberg/Blue M™
Forma®	Napco
Harris®	Precision™
Heraeus®	Puffer Hubbard®
Hybaid	Revco®
IEC®	Savant®
Jewett®	Sorvall
Jouan	

PRODUCT FAMILIES

Biosafety Cabinets
 Centrifuges
 CO₂ Incubators
 Cryopreservation
 Environmental Chambers
 Freeze Dryers
 Furnaces
 General Purpose Incubators
 Refrigerators and Freezers
 Ultra-low Temperature Freezers
 Microplate Instrumentation
 Ovens and Vacuum Ovens
 Shakers, Stirrers and Mixers
 Vacuum Concentrators
 Water Baths
 Water Purification

Catalogue Number for Single Installation	Catalogue Number for Additional Installations*	Equipment
ULS_INSTAL_1	ULS_INSTAL_1_ADD	<ul style="list-style-type: none"> • Table and Benchtop Centrifuges • Fridges and Freezers (including ULTs) • Microbiological Incubators and Benchtop Ovens • Microplate Washers and Microplate Incubators
ULS_INSTAL_2	ULS_INSTAL_2_ADD	<ul style="list-style-type: none"> • Floor Standing Ovens and Incubators • CO₂ Incubators • Benchtop and Floor Standing Centrifuges • Microplate Dispensers • Microplate Readers • Nucleic Acid Purification Instruments
ULS_INSTAL_3	ULS_INSTAL_3_ADD	<ul style="list-style-type: none"> • Type 1 Water Purification Systems
ULS_INSTAL_4	ULS_INSTAL_4_ADD	<ul style="list-style-type: none"> • Type 2 Water Purification Systems
ULS_INSTAL_5	ULS_INSTAL_5_ADD	<ul style="list-style-type: none"> • Clean Air Benches and Laminar Flow Cabinets
ULS_INSTAL_6	ULS_INSTAL_6_ADD	<ul style="list-style-type: none"> • Biosafety Cabinets Level 1
ULS_INSTAL_7	ULS_INSTAL_7_ADD	<ul style="list-style-type: none"> • Biosafety Cabinets Level 2 with Filter
ULS_COMMS_1	ULS_COMMS_1_ADD	<ul style="list-style-type: none"> • Commissioning of Biosafety Cabinets Level 1
ULS_COMMS_2	ULS_COMMS_2_ADD	<ul style="list-style-type: none"> • Commissioning Of Biosafety Cabinets Level 2 with Filter

* Additional installations must occur at same time as initial single installation

- Installation provided for other laboratory equipment, price on application
- Installation
 - Includes labor and travel
 - Connection to prepared service
 - Testing of equipment functionality and safety checks to manufacturers specification
 - Setup of user operating profiles on equipment
 - End-user training and general equipment maintenance
 - Single-point parameter verification (as appropriate for equipment)
- Commissioning for Biosafety Cabinets
 - Includes labor and travel
 - Full commissioning in accordance with BS EN 12469
 - Includes test of all HEPA filters for filter and filter seal integrity
 - Includes airflow check, alarm and electrical safety checks
 - Full validation certificate provided
 - Additional option to include test for operator protection (e.g., KI, DEHS according to local regulations)

ACT NOW!

Ask your sales representative to add installation service to your equipment quote.

Global Service Solutions

unitylabservices.com

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Please consult your local sales representative for details.

Austria
unity.at@thermofisher.com

Belgium
unity.bnl@thermofisher.com

Canada
unity.ca@thermofisher.com

Denmark
unity.dk@thermofisher.com

Finland
unity.fi@thermofisher.com

France
unity.fr@thermofisher.com

Germany
unity.de@thermofisher.com

Ireland
unity.ie@thermofisher.com

Italy
unity.it@thermofisher.com

Netherlands
unity.bnl@thermofisher.com

Norway
unity.no@thermofisher.com

Portugal
unity.pt@thermofisher.com

Spain
unity.es@thermofisher.com

Sweden
unity.se@thermofisher.com

Switzerland
unity.ch@thermofisher.com

UK
unity.uk@thermofisher.com

USA
unity.usa@thermofisher.com

Unity Lab Services

A Thermo Fisher Scientific Brand